

## MSH : Message Header

Location	Data Element	Data	Categorization
MSH-1	Field Separator		
MSH-2	Encoding Characters	^~\&	
MSH-3	Sending Application		
MSH-3.1	Namespace ID	NISTIISAPP	
MSH-3.2	Universal ID		
MSH-3.3	Universal ID Type		
MSH-4	Sending Facility		
MSH-4.1	Namespace ID	NISTIISFAC	
MSH-4.2	Universal ID		
MSH-4.3	Universal ID Type		
MSH-5	Receiving Application		
MSH-5.1	Namespace ID	NISTEHRAPP	
MSH-5.2	Universal ID		
MSH-5.3	Universal ID Type		
MSH-6	Receiving Facility		
MSH-6.1	Namespace ID	NISTEHRFAC	
MSH-6.2	Universal ID		
MSH-6.3	Universal ID Type		
MSH-7	Date/Time Of Message		
MSH-7.1	Time	20150625072816.602-0500	
MSH-9	Message Type		
MSH-9.1	Message Code	ACK	
MSH-9.2	Trigger Event	V04	
MSH-9.3	Message Structure	ACK	
MSH-10	Message Control ID	NIST-IZ-AD-10.2_Receive_ACK_Z23	
MSH-11	Processing ID		
MSH-11.1	Processing ID	P	
MSH-12	Version ID		
MSH-12.1	Version ID	2.5.1	
MSH-15	Accept Acknowledgment Type	NE	
MSH-16	Application Acknowledgment Type	NE	
MSH-21	Message Profile Identifier		
MSH-21.1	Entity Identifier	Z23	
MSH-21.2	Namespace ID	CDCPHINVS	
MSH-21.3	Universal ID		

MSH-21.4	Universal ID Type		
MSH-22	Sending Responsible Organization		
MSH-22.1	Organization Name	NISTIISFAC	
MSH-22.6	Assigning Authority		
MSH-22.6.1	Namespace ID		
MSH-22.6.2	Universal ID		
MSH-22.6.3	Universal ID Type		
MSH-22.7	Identifier Type Code		
MSH-22.10	Organization Identifier		
MSH-23	Receiving Responsible Organization		
MSH-23.1	Organization Name	NISTEHRFAC	
MSH-23.6	Assigning Authority		
MSH-23.6.1	Namespace ID		
MSH-23.6.2	Universal ID		
MSH-23.6.3	Universal ID Type		
MSH-23.7	Identifier Type Code		
MSH-23.10	Organization Identifier		

### MSA : Message Acknowledgment

Location	Data Element	Data	Categorization
MSA-1	Acknowledgment Code	AR	
MSA-2	Message Control ID	NIST-IZ-AD-10.1_Send_V04_Z22	

### ERR : Error

Location	Data Element	Data	Categorization
ERR-2	Error Location		
ERR-2.1	Segment ID	MSH	
ERR-2.2	Segment Sequence	12	
ERR-2.3	Field Position	1	
ERR-2.4	Field Repetition	1	
ERR-2.5	Component Number		
ERR-2.6	Sub-Component Number		
ERR-3	HL7 Error Code		
ERR-3.1	Identifier	203	
ERR-3.2	Text	Unsupported Version ID	

ERR-3.3	Name of Coding System	HL70357	
ERR-4	Severity	E	
ERR-5	Application Error Code		
ERR-5.1	Identifier		
ERR-5.2	Text		
ERR-5.3	Name of Coding System		
ERR-8	User Message	Version ID not recognized - message rejected	